



DRAWN BY:

DM

JUNKIN SAFETY APPL. CO.
3121 MILLERS LN
LOUISVILLE KY 40216

PT # **CT-150 CHANGING TABLE**

DATE 7/12/07

MATERIAL See Bill Materials

CT-150 ASSEMBLY CHART

Parts List			
ITEM	QTY	PART NUMBER	DESCRIPTION
1	4	CH-100-01	CORNER CASTING
2	2	CT-150-07	LOWER END TUBE
3	4	CT-150-08	UPRIGHT TUBE
4	4	MC-100-47	T-CASTING-POST
5	2	CT-150-10	UPPER CROSSBAR
6	2	CT-150-01	CAM
13	2	CT-150-14	KNOB
14	4	CT-150-09	INNER ELEVATING TUBE
15	2	CT-150-06	LOWER SIDE TUBE
17	8	CT-150-05	TABLETOP FRAME
24	4	JAS-800-13	SIDEARM
30	2	CT-150-12	Swivel Caster with Brake
31	1	CT-150-13	Swivel Caster Directional Lock
32	8	1/4 Flat Washer SAE	Washer B
33	8	MC-100-116	1/4-20 X 1-3/4 Hex Hd Cap Screw
34	8	CH-339	.250 WASHER
36	4	CT-150-25	1/4-20 X 1-1/2 Hex Cap Screw
37	4	CH-300-349	1/4-20 X 5/8 Long Socket hd Cap Screw
38	4	MC-100-131	1/4-20 NYLON LOCKING NUT
39	1	CT-150-03	TABLE TOP
40	2	CT-150-21	TABLETOP GUSSETT
41	2	CT-150-22	LOWER SIDE TUBE-OD