



**Note: PT# JSA-800-32L & 32R Swing Hook Spring
are Not Shown. It is assembled on item
"1" which is PT# JSA-800-24
(Swing Hook Spacer).**

DRAWN BY:

**JUNKIN SAFETY APPL. CO.
3121 MILLERS LN
LOUISVILLE KY 40216**

PT # **JSA-800-EH Pg 1**

DATE
6/22/11

MATERIAL

Parts List			
ITEM	QTY	PART NUMBER	DESCRIPTION
1	1	JSA-800-24	SWING HOOK SPACER
2	10	MC-100-103	Hex Nylon Locking Nut
3	1	JSA-800-34	ALL THREAD-LOWER CROSSBAR
4	1	JSA-800-20R	EXTENSION BRACKET
5	3	JSA-800-62	Regular Protruding Head Style Break Mandrel Blind Rivet - Inch
6	1	JSA-800-21R	SWING HOOK
7	1	JSA-800-19	Lower Extension Tube
8	1	JSA-800-60	SIDE TUBE STOP-UHMW
9	1	JSA-800-64	Push Button
10	1	JSA-800-18	SIDETUBE-EH
11	1	JSA-800-61	INNER TUBE GLIDE
12	2	JSA-800-65	10-24 x 3/8 Pan Hd Phillips
13	1	JSA-800-66	1/4-20 x 3.00 Hex Hd Bolt
14	1	JSA-800-28	ARM SPACER-LONG
15	1	JSA-800-67	1/4-20 x 1.50 Hex Bolt
18	2	JSA-800-68	1/4 INCH WASHER

DRAWN BY:

JUNKIN SAFETY APPL. CO.
3121 MILLERS LN
LOUISVILLE KY 40216

PT #

JSA-800-EH Pg 2

DATE

MATERIAL